

Chemblend® 3-Part Developer

***A World-Wide
Commitment
to Imaging
Excellence***

Chemblend® developer produces consistent quality on all types of film and in all makes of processors and cycles. A low fog level with balanced speed and contrast combine for *maximum detail*. Chemblend developer also has a 30% longer shelf life and runs clean in all processors, which means less maintenance and repairs.



CBM20
Chemical
Mixing
Station



Compare **Chemblend developer** to its competition—with decades of detail, it's tried and true.

Chemblend
Stock Number S60-040

Chemblend developer makes 5, 10, 15, and 20 gallons. Four bottles each part A, B, and C. For use with automatic processors.

Chemblend also comes in bulk form in an economical powder or liquid mix. Each makes 200 gallons.

Chemblend
Stock Number S60-047

Chemblend developer makes 5 and 10 gallons. Four bottles each part A, B, and C.



www.hrsimon.com

H.R. Simon and Company, Inc.
3515 Marmenco Court
Baltimore, Maryland 21230 • USA

TOLL FREE 800-638-9460
BALTIMORE 410-636-5555
FAX 410-636-4415
E-MAIL sales@hrsimon.com

S60-047 Chemblend 3 - Part Developer 10 Gallon Pack 32 lbs / 14.5 kg
S60-040 Chemblend 3 - Part Developer 20 Gallon Pack 61 lbs / 27.7 kg

TO MAKE		USE THIS AMOUNT OF WATER			USE THESE UNITS		
Gallons	Liters	Gallons +	Ounces	Liters	Part A	Part B	Part C
5	19	3	83	13.8	1	1	1
10	38	7	38	27.6	2	2	2
15	57	10	121	41.4	3	3	3
20	76	14	76	55.2	4	4	4

S50-090 Chemblend 3 - Part Developer 200 Gallon Bulk Liquid Mix 609 lbs / 277 kg

TO MAKE		USE THIS AMOUNT OF WATER			USE THESE GALLONS		
Gallons	Liters	Gallons +	Ounces	Liters	Part A	Part B	Part C
200	757	145	64	550.7	50	2.5	2

*Never leave containers half-full

Also Available: HS Automatic Developer 200 Gallon Bulk Powder Mix

PRODUCTS • Chemblend Developers • SEF² Developer • Autochem® Dental • HS Manual Solutions
• CBM20 Mixing Station • Imaging Services • Equipment and Accessories