

# Chemblend® Fixer

**3 Reasons to Buy  
Chemblend®  
Fixer...**  
**Runs Cleaner,  
Less Odor,  
Less Cost**



We pride ourselves in producing quality, providing our customers the best value at the best price.

The benefits of **Chemblend® Fixer** are:

- Chemblend runs cleaner – tests show it has the industry's lowest sulphide build-up characteristics.
- Highest archival standards
- Odorless formula makes the work environment more pleasant
- Economical

Your satisfaction with our products and our service is always our main focus.



CBM20  
Chemical  
Mixing  
Station

**Chemblend® Fixer**  
Stock Number S60-042

*Chemblend fixer for use with automatic processors. Each case makes up to 20 gallons.*

**Chemblend® Fixer**  
Stock Number S60-049

*Chemblend fixer for use with automatic processors. Each case makes up to 10 gallons.*

**Chemblend® Fixer**  
Stock Numbers S60-080  
and S60-081

*Chemblend automatic fixer makes 250 gallons. Packaged in one returnable 55 gallon drum and one re-usable 6 gallon DOT approved jug.*



www.hrsimon.com

**H.R. Simon and Company, Inc.**  
3515 Marmenco Court  
Baltimore, Maryland 21230 • USA

TOLL FREE 800-638-9460  
BALTIMORE 410-636-5555  
FAX 410-636-4415  
E-MAIL sales@hrsimon.com

S60-049		Chemblend Fixer 10 Gallon Pack			27 lbs / 12.3 kg	
S60-042		Chemblend Fixer 20 Gallon Pack			51 lbs / 23.2 kg	
TO MAKE		USE THIS AMOUNT OF WATER			USE THESE UNITS	
Gallons	Liters	Gallons +	Ounces	Liters	Part 1	Part 2
5	19	3	112	14.7	1	1
10	38	7	96	29.3	2	2
15	57	11	80	44	3	3
20	76	15	64	58.7	4	4
S60-080		Chemblend Fixer 250 Gallon Bulk Mix			609 lbs / 277 kg	
TO MAKE		USE THIS AMOUNT OF WATER			USE GAL / OUNCES	
Gallons	Liters	Gallons +	Ounces	Liters	Part 1	Part 2
50	189	38	96	146.7	10	1 / 32
100	378.6	77	64	293.3	20	2 / 64
150	568	116	32	440.1	30	3 / 96
200	757	155	0	586.8	40	5 / 0
250	946	193	96	733.5	50	6 / 32

Also Available: HS Automatic Fixer 300 Gallon Bulk Mix

PRODUCTS • Chemblend Developers • SEF<sup>2</sup> Developer • Autochem® Dental • HS Manual Solutions  
• CBM20 Mixing Station • Imaging Services • Equipment and Accessories